



More information on the website
radwag.com/ja/info,w1,DVE

5Y.10.PM 上皿天秤, 5Y.35.PM 上皿天秤, 5Y.60.PM 上皿天秤



5Y.10.PM 上皿天秤



5Y.35.PM 上皿天秤



5Y.60.PM 上皿天秤

The drawings, photos and graphics used are for illustrative purposes only.

機能



Autotest



Dosing



Percent Weighing



Parts counting



Peak hold



Formulation



Newton unit
measurement



Statistics



Checkweighing



IR sensors



GLP Procedures



Animal weighing



Pipettes Calibration



Air density correction



Density determination



Differential weighing



Ambient conditions
monitoring



Statistical Quality Control



Packaged Goods Control



ALIBI Memory



Wi-Fi



Under-pan weighing:
- 5Y.10.PM 上皿天秤
- 5Y.35.PM 上皿天秤

データシート

	5Y.10.PM 上皿天秤 WL-224-0001	5Y.35.PM 上皿天秤 WL-224-0005	5Y.60.PM 上皿天秤 WL-224-0007
計測パラメータ			
最大ひょう量 [Max]	10 kg	35 kg	60 kg
最小荷重	0,5 g	5 g	10 g
プリロード範囲	1 kg	3,5 kg	6 kg
最小表示 [d]	0,01 g	0,1 g	0,1 g
法定計量時の単位 [e]	0,1 g	1 g	1 g
風袋範囲	-10 kg	-35 kg	-60 kg
標準繰り返し性 [5% Max]	0,004 g	0,04 g	0,041 g
標準繰り返し性 [Max]	0,01 g	0,05 g	0,15 g
標準最小計量値 (USP)	8,2 g	100 g	82 g
標準最小計量値 (U=1%,k=2)	0,82 g	10 g	8,2 g
直線性	±0,025 g	±0,3 g	±0,3 g
安定化時間	1 s	1 s	1 s
調整	規格に準拠 (着脱式アダプター)	規格に準拠 (着脱式アダプター)	規格に準拠 (着脱式アダプター)
OIMLクラス	II	II	II
物理的パラメータ			
レベリングシステム	- LevelSENSING	- LevelSENSING	- LevelSENSING
ディスプレイ	10" メンブレン式 + バッテリー 赤外線式	10" メンブレン式 + バッテリー 赤外線式	10" メンブレン式 + バッテリー 赤外線式
納品構成	...	...	...
計量皿寸法	200×185 mm	350×260 mm	400×500 mm
梱包寸法	520×520×280 mm	520×520×280 mm	710×610×240 mm
正味重量	12,96 kg	12 kg	22 kg
総重量	14,2 kg	16 kg	26,5 kg
工事中			
保護クラス	IP 43	IP 43	IP 43
通信インターフェース			
通信インターフェース	2×USB-A, USB-C, RS 232 (COM3), HDMI, Ethernet, Wi-Fi, Hotspot	2×USB-A, USB-C, RS 232 (COM3), HDMI, Ethernet, Wi-Fi, Hotspot	2×USB-A, USB-C, RS 232 (COM3), HDMI, Ethernet, Wi-Fi, Hotspot
電氣的パラメータ			
電源	Adapter: 100 – 240V AC 50/60Hz 1A; 15V DC 2,4A : 12 – 15V DC 1,6A max; 10–19W*	Adapter: 100 – 240V AC 50/60Hz 1A; 15V DC 2,4A : 12 – 15V DC 1,6A max; 10–19W*	Adapter: 100 – 240V AC 50/60Hz 1A; 15V DC 2,4A : 12 – 15V DC 1,6A max; 10–19W*
環境条件			
動作温度	+10 – +40 °C	+10 – +40 °C	+10 – +40 °C

* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.



Additional fee for verification



アクセサリ (Additional Fee)

メディアボックス
RFIDタグ
防振台
電源アダプタ

底部吊り下げ計量
RS 232、RS 485 ケーブル
RS 232、RS 485 ケーブル
天秤用保護カバー

追加モジュール
プロ仕様の計量テーブル
天秤用保護カバー
バーコードスキャナ
天秤保管ケース

スタンド、壁取り付けキット、取り付けブラケット
レシートプリンタ
指紋リーダー
RS 232 – USB コンバータ

ソフトウェア (Additional Fee)

E2Rシステム
RADWAG Remote Desktop

Label Editor R02
Scales Editor 2.1

デバイスの寸法

5Y.10.PM 上皿天秤





